



Quick Western Transfer Buffer [10x] (w/o methanol)

Quick Western Transfer Buffer [1x] is an efficient, safe and non-toxic protein transfer product. The heat generated during the transfer process is less than that of ordinary transfer buffer, and no additional ice-water bath is needed. If the detection time is long, just use an ice box to cool down in the transfer tank. It can quickly complete the entire film transfer process within 5-30 minutes.

This product is safe and non-toxic and only needs to be diluted to 1x with water and absolute ethanol before use.

Catalog No.	503014
Size	100mL / 500mL
Product Category	Western-Blot
Storage/Stability	2 ~ 8°C/1 year
Shipping	Ambient