

## REAL-Detect Chlamydia Trachomatis qPCR Kit

REAL-Detect Chlamydia Trachomatis qPCR Kit uses qualitative fluorescence PCR method for invitro detection. The kit is supplied with reaction mix containing specific primers & fluorescent probes along with reaction enzyme & dNTPs; reverse transcriptase; positive & negative controls; and re-dissolved diluent. The kit is available in lyophilized format (stable at 25-30°C for 2 weeks).

Kit is compatible with all the standard Real-Time PCR systems, including:

ABI: 7000, 7300, 7500, 7900, Step One, Step One Plus

Agilent: Mx3000P, 3005P

Bio-Rad: CFX 96, iCycler iQ4, iCycleriQ5

Cepheid: Smart Cycler-II

MJ: Option2, Chromo4

Qiagen: Rotor Gene 6000

Roche: LightCycler680

Catalog No.	916192
Size	48 Test
Product Category	Bacterial Pathogen PCR
Storage/Stability	-20°C / 1 year
Shipping	Gel Packs

www.realgenelabs.com

For Research Use only